

5

TECHNICAL REPORT RK-83-4



LIQUID CHROMATOGRAPHIC ANALYSIS OF NITROCELLULOSE-BASE  
PROPELLANTS - FINAL REPORT OF FIRST YEAR EFFORT

Dr. James G. Carver  
Propulsion Directorate  
US Army Missile Laboratory

28 MARCH 1983

AD A129005



**U.S. ARMY MISSILE COMMAND**

*Redstone Arsenal, Alabama 35809*

*Approved for public release; distribution unlimited.*

DTIC FILE COPY

DTIC  
ELECTE  
JUN 7 1983  
S A D

#### **DISPOSITION INSTRUCTIONS**

**DESTROY THIS REPORT WHEN IT IS NO LONGER NEEDED. DO NOT  
RETURN IT TO THE ORIGINATOR.**

#### **DISCLAIMER**

**THE FINDINGS IN THIS REPORT ARE NOT TO BE CONSTRUED AS AN  
OFFICIAL DEPARTMENT OF THE ARMY POSITION UNLESS SO DESIGN-  
ATED BY OTHER AUTHORIZED DOCUMENTS.**

#### **TRADE NAMES**

**USE OF TRADE NAMES OR MANUFACTURERS IN THIS REPORT DOES  
NOT CONSTITUTE AN OFFICIAL INDORSEMENT OR APPROVAL OF  
THE USE OF SUCH COMMERCIAL HARDWARE OR SOFTWARE.**

UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

REPORT DOCUMENTATION PAGE		READ INSTRUCTIONS BEFORE COMPLETING FORM	
1. REPORT NUMBER RK-83-4	2. GOVT ACCESSION NO. <b>A129003</b>	3. RECIPIENT'S CATALOG NUMBER	
4. TITLE (and Subtitle) Liquid Chromatographic Analysis of Nitrocellulose-Base Propellants - Final Report of First Year Effort		5. TYPE OF REPORT & PERIOD COVERED Technical Report	
7. AUTHOR(s) Dr. James G. Carver		6. PERFORMING ORG. REPORT NUMBER	
9. PERFORMING ORGANIZATION NAME AND ADDRESS Commander, US Army Missile Command ATTN: DRSMI-RK Redstone Arsenal, AL 35898		8. CONTRACT OR GRANT NUMBER(s)	
11. CONTROLLING OFFICE NAME AND ADDRESS Commander, US Army Missile Command ATTN: DRSMI-RPT Redstone Arsenal, AL 35898		10. PROGRAM ELEMENT, PROJECT, TASK AREA & WORK UNIT NUMBERS	
14. MONITORING AGENCY NAME & ADDRESS (if different from Controlling Office)		12. REPORT DATE 28 March 1983	
		13. NUMBER OF PAGES 29	
		15. SECURITY CLASS. (of this report) UNCLASSIFIED	
		15a. DECLASSIFICATION/DOWNGRADING SCHEDULE	
16. DISTRIBUTION STATEMENT (of this Report) Approved for public release; distribution unlimited.			
17. DISTRIBUTION STATEMENT (of the abstract entered in Block 20, if different from Report)			
13. SUPPLEMENTARY NOTES			
19. KEY WORDS (Continue on reverse side if necessary and identify by block number) Liquid Chromatography                      Method Development Propellant Analysis Nitrocellulose Base Propellants			
20. ABSTRACT (Continue on reverse side if necessary and identify by block number) A procedure has been developed that is able to predict the optimum conditions for separation of the ingredients of a complex mixture by liquid chromatography. The procedure involves establishing a data base of retention times for a large number of propellant ingredients in 7 predetermined combinations of 3 solvents. When a new combination of some of these ingredients is encountered, a computer analysis of the retention times of the ingredients from the data base will predict the best solvent composition to separate all the ingredients. The details of the procedure and several examples of its use are discussed.			

DD FORM 1473

JAN 73

EDITION OF 1 NOV 65 IS OBSOLETE

UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE(When Data Entered)

ABSTRACT - Continued

This project has been accomplished as part of the US Army Materials Testing Technology Program, which has for its objective the timely establishment of testing techniques, procedures, or prototype equipment (in mechanical, chemical or nondestructive testing) to insure efficient inspection methods for materiel/ material procured or maintained by DARCOM.

SECURITY CLASSIFICATION OF THIS PAGE(When Data Entered)

# TABLE OF CONTENTS

	<u>Page No.</u>
I. INTRODUCTION .....	3
II. EXPERIMENTAL .....	3
A. EQUIPMENT .....	3
B. SAMPLE PREPARATION .....	4
III. PROCEDURE .....	4
IV. RESULTS .....	6
V. CONCLUSION .....	7
DISTRIBUTION .....	20

DISTRIBUTION FOR	
ALL	<input checked="" type="checkbox"/>
SECRET	<input type="checkbox"/>
CONFIDENTIAL	<input type="checkbox"/>
DISTRIBUTION	
Availability Codes	
AVAIL and/or	Special
A	

010 100 010

## LIST OF ILLUSTRATIONS

	<u>Page Number</u>
Figure 1. Interaction of Solvents .....	10
Figure 2. Response Surface for Gunpowder Ingredients .....	11
Figure 3. Chromatogram of Mixture of Ingredients .....	12
Figure 4. Gunpowder A Analysis .....	13
Figure 5. Gunpowder B Analysis .....	14
Figure 6. Gunpowder C Analysis .....	15
Figure 7. Response Surface for Propellant Ingredients .....	16
Figure 8. Chromatogram of Propellant .....	17
Figure 9. Chromatogram of Propellant .....	18
Figure 10. Chromatogram of Propellant .....	19

## I. INTRODUCTION

Nitrocellulose-based propellants are widely used throughout the military and civilian areas in guns and rockets. There are almost as many different procedures for analysis as there are propellants. Even in the area of liquid chromatography, frequently, more than one procedure is used in the analysis of the propellant. Most techniques have been designed to answer a specific question about the purity or amount of two or three ingredients. However, most propellants have at least three or four original ingredients and frequently one or two other reaction products that can be easily extracted.

Since these ingredients can be extracted into a solvent it would seem logical to try to analyse for all them at once. A major difficulty to this in the past has been that silica columns are difficult to work with. Precise control of water, alcohols, and other polar ingredients is necessary to obtain reproducible retention times and amounts. Reverse phase columns-octadecylsilane or C-18 do not have this problem but investigators have expressed difficulty in obtaining good separations of more than a few ingredients with the binary mixtures commonly used.

Recently some significant interest has been observed in using more than 2 solvent systems with C-18 columns. Manufacturers have started producing HPLC equipment that can pump 3 or 4 solvents simultaneously. The reason for a multi-solvent system is that separations are accomplished by more than differences in polarity. Snyder<sup>1</sup> has proposed that solvents should be classified as to three properties - proton acceptor, proton donor, and dipole interaction as indicated in Figure 1. It is believed that these properties make the major contributions to separations in HPLC. Kirkland<sup>2</sup> has developed a procedure consisting of 7 experiments which can be used to predict the best solvent composition for separating ingredients in a mixture. A program was undertaken to develop a procedure that will separate the maximum number of commonly used propellant ingredients possible.

## II. EXPERIMENTAL

### A. Equipment

All solvents were obtained from Burdick and Jackson Laboratories. Prior to use the solvents were filtered and degassed in an ultrasonic bath. Water was distilled and deionized by a MILLI-Q system from Millipore. The liquid chromatographic system consisted of three M6000A pumps, a 720 system controller, a 710B WISP sample injector, and RCM 100 radial compression system, all from Water Associates. The column used was a Waters Associates C-18 radial PAC column using 10 micron silica coated with octadecyl silane. The detector was a LC 85 variable wavelength UV detector from Perkin-Elmer. Responses were most often monitored at 215 nm. The signal from the detector was analyzed by a Hewlett-Packard 3353E laboratory automation system.

---

<sup>1</sup> Snyder, L. R. J Chrom Sci 16, (1978), 223.

<sup>2</sup> Glatch, J. L., Kirkland, J. J., Squire, K. M., and Miner, J. M.; J. Chromatography, 199, (1980), 57.

## B. Sample Preparation

The missile propellant was cut into small pieces and approximately 400 mg were accurately weighed out and placed in a test tube. To this was added an accurately weighed amount of diethylphthalate as an internal standard. Twenty-five ml of acetonitrile was added to the test tube. After one hour the contents of the test tube were subjected to ultrasonic vibrations using a Sonicator micro tip for 15 minutes. The solution was then filtered and placed in a WISP sample vial for analysis.

The gun powder was prepared in a similar manner except that the acetonitrile was diluted and there was no soak time. The reason from this modification was that the powder was not cross linked and tended to form a gum in pure acetonitrile. Usually 50% acetonitrile in water was sufficient to effect complete extraction.

## III. PROCEDURE

There are numerous papers and books on ways to evaluate mixtures in general. The paper by Snee<sup>3</sup> presents a good summary of many of the techniques in a manner easily understood by the uninitiated. Examination of these methods in light of Snyder's<sup>1</sup> paper permits development of an experimental procedure involving only 7 to 10 experiments to obtain a complete picture of the separation properties of the three solvents and their mixtures. Snyder<sup>1</sup> suggested the use of methanol, acetonitrile, and tetrahydrofuran as the three solvents of choice for reversed phase chromatography. In the real world we want a separation or resolution of the ingredients in a sample within a certain time limit. It has been shown<sup>4</sup> that resolution will increase with

$\frac{K}{1+K}$  but increases in values of K greater than 10 or 15 have only minimal influence on the resolution. The first step in developing a separation is to adjust the three solvents strengths by dilution with an inert carrier, such as water, until the retention times of all the ingredients in the sample are between 2 and 30 minutes (assuming non retained peaks elute at 1 minute). This produces 3 solvents of approximately the same strengths but greatly differing selectivities. Once the dilution of one solvent has been selected, the desired dilution of the other two can be roughly estimated as described by Kirkland<sup>5,2</sup>. For methanol, acetonitrile and tetrahydrofuran the following calculations were made:

$$\phi_a P_a + \phi_b P_b = \phi_A P_A + \phi_B P_B \quad (1)$$

Where  $\phi_A$  = % by vol of A

<sup>3</sup> Snee, R. D.; Chemtech; Nov (1979) 702.

<sup>4</sup> "Introduction to Modern Liquid Chromatography", L. R. Snyder and J. J. Kirkland, Wiley, New York, 2nd edition, 1979, Ch 2.

<sup>5</sup> Reference 4, Chapter 6.



$$P_{H_2O} = 10.2 \quad \text{for ACN: } 50\%(5.8) + 50\%(10.2) = 8.00$$

$$P_{ACN} = 5.8 \quad \text{for THF: } 38\%(4.0) + 64\%(10.2) = 7.97$$

$$P_{MEOH} = 5.1 \quad \text{for MECH: } 43\%(5.1) + 57\%(10.2) = 8.01$$

$$P_{THF} = 4.0$$

The dilutions finally selected were 60% methanol in water, 50% acetonitrile in water, and 42% tetrahydrofuran in water. For convenience these will be identified as MECH, ACN, and THF.

At this point the 7 to 10 experiments suggested by Snee can be performed. Here the retention times of each ingredient are determined in each of the mixtures of the three solvents as shown in Table I. Only the first 7 are necessary to generate the response surface but mixtures 8, 9, and 10 can be used to check the accuracy of the analysis.

TABLE I. SOLVENT COMPOSITIONS

Experiment	1	2	3	4	5	6	7	8	9	10
MECH	100	0	0	50	50	0	33.3	16.6	16.6	66.4
ACN	0	100	0	50	0	50	33.3	16.6	66.4	16.6
THF	0	0	100	0	50	50	33.3	66.4	16.6	16.6

Once the retention times have been obtained the data can be analyzed to develop a response surface for the separations. The resolution of any two ingredients can be calculated for any solvent mix of composition ( $X_1A + X_2B + X_3C$ ) from equation (2)

$$R = B_1X_1 + B_2X_2 + B_3X_3 + B_4X_1X_2 + B_5X_1X_3 + B_6X_2X_3 + B_7X_1X_2X_3 \quad (2)$$

The values for  $B_1$  through  $B_7$  are derived from equations 3 through 9.

Here  $Y_1$  is the resolution of a peak pair in solvent mixture 1.

$$B_1 = Y_1 \quad (3)$$

$$B_2 = Y_2 \quad (4)$$

$$B_3 = Y_3 \quad (5)$$

$$B_4 = 4Y_4 - 2(Y_1 + Y_2) \quad (6)$$

$$B_5 = 4Y_5 - 2(Y_1 + Y_3) \quad (7)$$

$$B_6 = 4Y_6 - 2(Y_2 + Y_3) \quad (8)$$

$$B_7 = 27*Y_7 - 12(Y_4 + Y_5 + Y_6) + 3(Y_1 + Y_2 + Y_3) \quad (9)$$

A computer program is being developed to perform these calculations and plot the response surface for selected resolutions of pairs of ingredients. The analysis is complicated by changes in order of elution of several ingredients as the solvent composition is varied. For this reason each ingredient in the mixture must be evaluated with each of the other ingredients for resolution.

#### IV. RESULTS

Samples of three gunpowders were obtained from a commercial manufacturer. A mixture was prepared that contained all of the ingredients and the probable reaction products. Table II lists the retention times of each ingredient for the 7 solvent mixtures suggested.

The B values calculated for several pairs of ingredients are listed in Table III.

Most of the ingredients are always resolved by a minimum value of 1.25 and are not considered critical to the separation. Figure 2 is a superimposed graph of the lines of resolution of the critical ingredients.

As can be seen in Figure 2 there is a small area of the graph where all ingredients can be resolved by the minimum value of 1.25. The solvent mixture in the middle of this area was selected and the 10 ingredient mixture was analysed. The resulting chromatograph, Figure 3, shows that indeed all the ingredients are well separated. At this point the three gun powder samples were prepared for analysis. Samples A and B did not contain Diphenylphthalate so this was selected as the internal standard. The internal standard for sample C was Dimethylphthalate since it already contained Diphenylphthalate. Chromatograms and averaged analysis of each sample are given in Figures 4, 5, and 6 and Tables IV, V, and VI.

A sample of cross linked double-base propellant was obtained for analysis. The individual ingredients were placed in 3 different mixtures to aid in identification and analysed by the same solvent scheme in Table I. The B coefficients for the critical peak pairs are given in Table VII and the superimposed response surfaces are given in Figure 7. Here a much larger area of solvent choice occurred. The propellant sample was prepared using Diethylphthalate as an internal standard. It was analysed using three different solvent mixtures within the indicated area and these are shown in Figures 8, 9, and 10.

The results of the analysis were in good agreement with other analyses of the propellant.

It is interesting to note that gas chromatography has difficulty in analysing for thermally labile ingredients, such as nitroglycerine, and low volatility materials, such as HMX. Liquid chromatography is generally limited to solubility and detectability of the ingredient. Nitroglycerine does not absorb well in the UV at 254 nm. To detect it one must go down to about 214 or 220 nm. HMX has a low solubility in most solvents. In the first sample

preparation of the propellant sample the HMX analysis was low while all other ingredients were accurate. Examination of the procedure disclosed that not enough solvent had been used to dissolve all the HMX. Subsequent preparation used larger amounts of solvent and the analysis improved.

## V. CONCLUSIONS

The procedure described here seems to have almost universal application in liquid chromatography. Once the retention times for a list of ingredients is developed no further experimentation is required. When a new propellant composition is encountered, the solvent composition necessary for optimum separation of the ingredients can be quickly predicted with the aid of a computer.

A data base of 30 of the most commonly used propellant ingredients is being prepared for further testing. Currently a mixture of 10 ingredients can be separated with one solvent mixture. It is the ultimate goal that all 30 ingredients can be separated with one solvent mixture but it is doubtful that this can be achieved.

A variable that has not been examined yet is the influence of using different columns. During the next year this will be evaluated by using two C18 columns from the same manufacturer and other C18 columns from different manufacturers. It is expected that the within manufacturer variations will be small and not require changes in the optimum solvent mixture. There will probably be significant differences between manufacturers but the extent of the differences can not be predicted at this time.

TABLE II. GUNPOWDER INGREDIENT RETENTION TIMES

SOLVENT MIX INGREDIENTS	1	2	3	4	5	6	7
DMP	1.70	1.86	1.24	1.89	1.33	1.33	1.42
DEP	3.09	3.99	2.01	3.94	2.23	2.33	2.57
2, 4 DNT	3.35	3.34	1.91	3.15	2.15	2.66	2.90
NG	2.70	2.03	3.01	2.53	3.13	5.42	3.52
EC	5.61	16.5	3.70	13.20	4.10	9.40	6.28
NDPA	4.84	6.24	5.67	6.25	2.74	5.23	4.94
DPA	5.61	6.90	5.67	7.43	5.69	5.79	7.89
DPP	13.96	21.0	5.26	22.5	7.13	10.86	11.74
DBP	6.99	35	7.92	32	10.15	15	16

DMP = DIMETHYL PHTHALATE  
 DEP = DIETHYL PHTHALATE  
 2,4 DNT = 2,4 DINITRO TOLUENE  
 NG = NITROGLYCERINE  
 EC = ETHYL CENTRALITE

NDPA = N-NITROSODIPHENYLAMINE  
 DPA = DIPHENYLAMINE  
 DPP = DIPHENYLPHTHALATE  
 DBP = DIBUTYLPHTHALATE

TABLE III. GUNPOWDER B COEFFICIENTS FOR PAIRS OF INGREDIENTS

	B1	B2	B3	B4	B5	B7	B7
DMP-NG	6	1.02	10.62	1.32	9.96	74.88	-77.04
DEP-2,4 DNT	1.56	-3.9	-.6	-14.28	-3.84	16.92	83.52
2,4 DNT-NG	-3.9	-7.86	6.6	8.64	18.12	68.76	-139.68
NG-NDPA	12.84	25.26	15.96	13.08	-66.96	-87	166.14
EC-NDPA	-4.62	-61.56	11.82	-34.44	-47.04	-.6	518.4
EC-DPA	0	-57.6	11.82	-23.28	14.52	4.92	684.36
NDPA-DPA	4.62	3.96	0	11.16	61.56	5.52	165.96
DPA-DPP	50.1	84.6	-2.46	92.28	-60.72	-42.6	-533.34

TABLE IV. GUNPOWDER A ANALYSIS

Nitroglycerine	$11.421 \pm .057$
DNT	$.354 \pm .001$
NDPA	$.327 \pm .009$
DPA	$.652 \pm .008$
2-NDPA	$.069 \pm .001$
DBP	$1.980 \pm .002$

TABLE V. GUNPOWDER B ANALYSIS

Nitroglycerine	$46.4 \pm .7$
DNT	$.234 \pm .002$
NDPA	$.294 \pm .003$
DPA	$.054 \pm .001$
2-NDPA	$.050 \pm .000$
DBP	$.115 \pm .004$

TABLE VI. GUNPOWDER C ANALYSIS

Nitroglycerine	$6.649 \pm .007$
DNT	$.067 \pm .006$
NDPA	$.264 \pm .003$
DPA	$.243 \pm .004$
EC	$.535 \pm .001$
2-NDPA	$.098 \pm .001$
DPP	$.592 \pm .001$
DEP	$.043 \pm .002$

TABLE VII. B COEFFICIENTS FOR PROPELLANT FOR PAIRS OF INGREDIENTS

	B1	B2	B3	B4	B5	B6	B7
1-2 BTTN-NG	-3.24	-5.54	-12.72	-5.76	24.48	9.6	-20.7
1-6 BTTN-4-NDPA	-6.42	-3.12	-15.36	16.68	9.24	-55.92	-26.1
1-8 BTTN-DEP	-.9	6.36	-18.72	.12	10.2	-76.08	34.74
2-4 NG-N-MNA	-3.18	2.28	-2.64	22.44	-15.24	-65.52	-5.4
2-8 NG-DEP	2.34	11.76	-6.0	5.88	-14.28	-85.68	55.44
3-4 MNA-N-MNA	2.1	1.92	5.64	21.72	-7.56	-2.88	-85.14
3-8 MNA-DEP	7.62	11.4	2.28	5.16	-6.6	-23.04	-24.3
4-8 N-MNA-DEP	5.52	9.48	3.36	16.56	.96	-20.16	60.84




Figure 1. Interaction of solvents.




Figure 2. Response surface for gunpowder ingredients.




Figure 3. Chromatogram of mixture of ingredients.






Figure 4. Gunpowder A analysis.




Figure 5. Gunpowder B analysis.




Figure 6. Gunpowder C analysis.




Figure 7. Response surface for propellant ingredients.




Figure 8. Chromatogram of propellant.




Figure 9. Chromatogram of propellant.




Figure 10. Chromatogram of propellant.

# DISTRIBUTION LIST

	<u>No. of Copies</u>
Commander US Army Materiel Development and Readiness Command ATTN: DRCQA, Mr. S. Lorber 5001 Eisenhower Avenue Alexandria, VA 22333	1
Commander US Army Materiel Development and Readiness Command ATTN: DRCQA, Mr. J. Lockerd 5001 Eisenhower Avenue Alexandria, VA 22333	1
Commander US Army Materiel Development and Readiness Command ATTN: DRCQA-EQ, Mr. K. Snelly 5001 Eisenhower Avenue Alexandria, VA 22333	1
Commander US Army Materiel Development and Readiness Command ATTN: DRCQA-EQ, Mr. H. Light 5001 Eisenhower Avenue Alexandria, VA 22333	1
Commander US Army Materiel Development and Readiness Command ATTN: DRCQA-P, Mr. R. Tiner 5001 Eisenhower Avenue Alexandria, VA 22333	1
Commander US Army Materiel Development and Readiness Command ATTN: DRCMT, Mr. F. Michel 5001 Eisenhower Avenue Alexandria, VA 22333	1
Commander US Army Materiel Development and Readiness Command ATTN: DRCMT, Mr. J. L. Baer 5001 Eisenhower Avenue Alexandria, VA 22333	1
Commander US Army Materiel Development and Readiness Command ATTN: DRCTM-S, Mr. T. Wolf 5001 Eisenhower Avenue Alexandria, VA 22333	1



	<u>No. of Copies</u>
Commander US Army Armament Materiel Readiness Command ATTN: DRSAR-QA, Mr. L. Griffin Rock Island, IL 61299	1
Commander U. S. Army Armament Materiel Readiness Command ATTN: DRSAR-QAE, Mr. R. Fer Rock Island, IL 61299	1
Commander U. S. Army Armament Materiel Readiness Command ATTN: DRSAR-QAS, Mr. J. W. Hawkins Rock Island, IL 61299	1
Commander U.S. Army Armament R&D Command ATTN: DRDAR-QA, Mr. H. Lazar Dover, NJ 07801	1
Commander U. S. Army Armament R&D Command ATTN: DRDAR-QAS-D, Mr. J. Vega Dover, NJ 07801	1
Commander U. S. Army Aviation R&D Command ATTN: DRDAV-Q, Mr. E. J. Hollman St. Louis, MO 63120	1
Commander U. S. Army Aviation R&D Command ATTN: DRDAV-QE, Mr. J. T. Conroy St. Louis, MO 63120	1
Commander U.S. Army Aviation R&D Command ATTN: DRDAV-QE, Mr. W. Baker St. Louis, MO 63120	1
Commander U. S. Army Depot Systems Command ATTN: DRSDE-Q, Mr. W. D. Henry Letterkenny Army Depot Chambersburg, PA 17201	1
Commander U. S. Army Communications-Electronics Command ATTN: DRSEL-PA, Mr. A. D'Angelo Fort Monmouth, NJ 07703	1

	<u>No. of Copies</u>
Commander U. S. Army Communications-Electronics Command ATTN: DRSEL-PA-M, Mr. C. Faulkner Fort Monmouth, NJ 07703	1
Commander U. S. Army Avionics R&D Activity ATTN: DAVAA-N, Mr. Matthew J. Fisher Fort Monmouth, NJ 07703	1
Commander U. S. Army Electronics R&D Command ATTN: DRDEL-PA, Mr. S. Alster 2300 Powder Mill Road Adelphi, MD 20783	1
Commander U. S. Army Electronics R&D Command ATTN: DRDEL-PA, Mr. J. Goon 2500 Powder Mill Road Adelphi, MD 20783	1
Commander Harry Diamond Laboratories ATTN: DELHD-PO-P, Mr. J. Hoke 2300 Powder Mill Road Adelphi, MD 20783	1
Commander Harry Diamond Laboratories ATTN: DELHD-DE-OP, Mr. A. B. Goldberg Adelphi, MD 20783	1
Director Night Vision & Electro-Optics Laboratory ATTN: DELNV-PA, Mr. D. Jenkins Fort Belvoir, VA 22060	1
Commander U. S. Army Missile Command ATTN: DRSMI-QET, Mr. D. Mcvey Redstone Arsenal, AL 35898	1
Commander U. S. Army Missile Command ATTN: DRSMI-RKP, Mr. J. Wright Redstone Arsenal, AL 35898	1

	<u>No. of Copies</u>
Commander U. S. Army Missile Command ATTN: DRSMI-RCK, Mr. B. J. Alley Redstone Arsenal, AL 35898	1
Commander U. S. Army Missile Command ATTN: DRSMI-RST, Mr. L. Chapman Redstone Arsenal, AL 35898	1
Commander U. S. Army Mobility Equipment Research and Development Command ATTN: DRDME-T, Mr. E. D. Hardin Fort Belvoir, VA 22060	1
Commander U. S. Army Mobility Equipment Research and Development Command ATTN: DRDME-T, Mr. R. Bolka Fort Belvoir, VA 22060	1
Commander U. S. Army Mobility Equipment R&D Command ATTN: DRDME-V, Mr. E. York Fort Belvoir, VA 22060	1
Commander U. S. Army Mobility Equipment R&D Command ATTN: DRDME-V, Mrs. Wells Fort Belvoir, VA 22060	1
U. S. Army Materiel System Analysis Activity ATTN: DRXSY-MP Aberdeen Proving Ground, MD 21005	1
Commander U. S. Army Missile Command ATTN: DRSMI-Q, Mr. J. E. Shepard Redstone Arsenal, AL 35898	1
Commander U. S. Army Mobility Equipment R&D Command ATTN: DRDME-VC, Mr. D. Emeric Fort Belvoir, VA 22060	1

	<u>No. of Copies</u>
Commander U. S. Army Armament R&D Command ATTN: DRDAR-QAR-I, Mr. M. H. Weinberg Dover, NJ 07801	1
Commander U. S. Army Armament R&D Command ATTN: DRDAR-SCM-P, Mr. J. Beetle Dover, NJ 07801	1
Commander Watervliet Arsenal ATTN: SARVV-QA, Mr. J. R. Germain Watervliet, NY 12189	1
Commander Watervliet Arsenal ATTN: SARVV-QA, Mr. J. Miller Watervliet, NY 12189	1
Commander Watervliet Arsenal ATTN: SARVV-QAE-S, Mr. S. Krupski Watervliet, NY 12189	1
Commander Watervliet Arsenal ATTN: SARVV-PPI, Mr. R. McCabe Watervliet, NY 12189	1
Commander U. S. Army Armament R&D Command ATTN: DRDAR-QAS-A, Mr. Geza Pap Dover, NJ 07801	1
Commander ARRADCOM, Edgewood Area ATTN: DRDAR-QAC-E, Dr. W. J. Maurits Aberdeen Proving Ground, MD 21005	1
Commander U. S. Army Natick R&D Command ATTN: DRDNA-E, Mr. E. F. Levell Kansas Street Natick, MA 01760	1

	<u>No. of Copies</u>
Commander U. S. Army Natick R&D Command ATTN: DRDNA-EML, Mr. F. Civilikas Kansas Street Natick, MA 01760	1
Commander U. S. Army Tank-Automotive Command ATTN: DRSTA-Q, Mr. L. Barnett Warren, MI 48090	1
Commander U. S. Army Tank-Automotive Command ATTN: DRSTA-QAT, Mr. D. Gamache Warren, MI 48090	1
Commander U. S. Army Tank-Automotive Command ATTN: DRSTA-QAT, Mr. N. Keck Warren, MI 48090	1
Commander U. S. Army Tank-Automotive Command ATTN: DRSTA-QAT, Mr. C. Kedzior Warren, MI 48090	1
Commander U. S. Army Tank-Automotive Command ATTN: DRSTA-RCK, Dr. J. Chevalier Warren, MI 48090	1
Commander U. S. Army Tank-Automotive Command ATTN: DRSTA-RCK, Mr. S. Goodman Warren, MI 48090	1
Commander U. S. Army Tank-Automotive Command ATTN: DRSTA-RCKM, Mr. D. Ostberg Warren, MI 48090	1
Commander U. S. Army Test & Evaluation Command ATTN: DRSTE-AD-M, R. G. Sheldon Aberdeen Proving Ground, MD 21005	1
Commander U. S. Army Test & Evaluation Command ATTN: DRSTE-AD-M, Mr. R. Huddleston Aberdeen Proving Ground, MD 21005	1

	<u>No. of Copies</u>
Commander Pine Bluff Arsenal ATTN: SARPB-QA, Mr. T. L. Ashcroft Pine Bluff, AR 71611	1
Commander Rocky Mountain Arsenal ATTN: SARRM-QA, Mr. S. G. Jacobs Commerce City, CO 80022	1
Commander Rock Island Arsenal ATTN: SARRI-QA, Mr. W. J. Betts, Jr. Rock Island, IL 61229	1
Commander Aberdeen Proving Ground ATTN: STEAP-G, Mr. C. Klarich Aberdeen Proving Ground, MD 21005	1
Commander U. S. Army Electronic Proving Ground ATTN: STEEP-MT, LTC Melvin Fowler Fort Huachuca, AZ 85613	1
Commander U. S. Army Dugway Proving Ground ATTN: STEDP-PO-Q Dugway, UT 84022	1
Commander U. S. Army Jefferson Proving Ground ATTN: STEJP-TD, Mr. I. K. Peterson Madelson, IN 47250	1
Commander U. S. Army Yuma Proving Ground ATTN: STEYP-MTD, LTC R. Gill'gar Yuma, AZ 85364	1
Commander U. S. Army White Sands Missile Range ATTN: STEWS-QA, Mr. G. F. Bigelow White Sands Missile Range, MN 88002	1
Commander US Army Tropic Test Center ATTN: STETC APO Miami, FL 34004	1

No. of  
Copies

Commander  
U. S. Army Troop Support and Aviation Materiel  
Readiness Command  
ATTN: DRSTS-Q, Mr. W. G. Creel  
4300 Goodfellow Boulevard  
St. Louis, MO 63120

1

Commander  
Anniston Army Depot  
ATTN: SDSAN-DQA, Mr. Pennington  
Anniston, AL 36202

1

Commander  
Corpus Christi Army Depot  
ATTN: SDSCC-Q, Mr. L. E. Ervin  
Corpus Christi, TX 78419

1

Commander  
Letterkenny Army Depot  
ATTN: SDSLE-Q, Mr. W. J. Walker  
Chambersburg, PA 17201

1

Commander  
New Cumberland Army Depot  
ATTN: SDSNC-Q, Mr. A. T. Holderbach  
New Cumberland, PA 17070

1

Commander  
Red River Army Depot  
ATTN: SDSRR-Q, Mr. W. D. Wertz  
Texarkana, TX 75501

1

Commander  
Lexington-Blue Grass Activity  
ATTN: SDSAX-LQ, Mr. P. Jackson  
Lexington, KY 40511

1

Commander  
Sacramento Army Depot  
ATTN: SDSSA-Q, Mr. C. H. Clawson  
Sacramento, CA 95801

1

Commander  
Savanna Army Depot  
ATTN: SDSSV-VE, LTC P. E. James  
Savanna, IL 61074

1

	<u>No. of Copies</u>
Commander Seneca Army Depot ATTN: SDSSE-R, Mr. P. W. Chavez Romulus, NY 14541	1
Commander Sharpe Army Depot ATTN: SDSSE-Q, Mr. J. E. Seyfried Lathrop, CA 95330	1
Commander Sierra Army Depot ATTN: SDSSI-QA, Mr. L. B. Lampe Barlong, CA 96113	1
Commander Tooele Army Depot ATTN: SDSTE-QA, Mr. R. M. Rich Tooele, UT 84074	1
Commander Pueblo Depot Activity ATTN: SDSTE-PUG, Mr. P. S. Coen Pueblo, CO 81001	1
Commander Defense Technical Information Center Cameron Station, Bldg # 5 5010 Duke St. Alexandria, VA 22314	5
Commander USA foreign Science & Tech Center ATTN: DRX31-SD3 220 7th St., N.E. Charlottesville, VA 22901	1
Office of Deputy Chief of Staff for Research, Development & Acq. ATTN: DAMA-ARZ-E DAMA-CSS Washington, D.C. 20310	1 1
Commander Army Research Office ATTN: Dr. Mayer Mr. J. J. Murray P.O. Box 12211 Research Triangle Park, NC 27709	1 1



	<u>No. of Copies</u>
Chemical Propulsion Information Agency John Hopkins University John Hopkins Rd Laurel, MD 20810	1
Commander Tobyhanna Army Depot ATTN: SDSTO-Q, Mr. W. J. Lord Tobyhanna, PA 18466	1
Director U. S. Army Industrial Base Eng. Activity ATTN: DRXIB-MT, Mr. D. Brin Rock Island, IL 61299	1
Director Army Materials and Mechanics Research Center ATTN: DRXMR-STQ	5
DRXMR-PL	2
Watertown, MA 02172	
Commander U. S. Army Missile Command ATTN: DRSMI-RPT	1
DRSMI-LP	1
DRSMI-RPR	15
Redstone Arsenal, AL 35898	
Commander US Army Materiel Systems Analysis Activity ATTN: DRXSY-MP Aberdeen Proving Ground, MD 21005	1